Clean up with Cingular's most cost effective Pulse Lavage on the market.

$\mathbf{C} \sim \mathbb{Pulse}^{\mathsf{m}}$

ALL YOU NEED IS RIGHT HERE

Introducing Cingular's C~PulseTM Disposable Pulse Lavage.

Cingular's new generation C~PulseTM, with FIVE unique features make it user to use, eliminating user risks with trigger position and battery housing. Independent tubing and flow regulator maximising the user experience – guaranteed, sterile, saving hospitals valuable and expensive operating room time and costs - the Cingular Pledge.


INDICATION OF USE

- Removing blood and debris.
- Supplying an improved bone-to-cement interface.
- Reducing the risk of fat embolism.

SPECIFICATION


POWER: 8 x 1.5V CELL

FLOW RATE: 1200ml/min HIGH 600ml/min LOW



C - Pulse

COMPATIBLE TIPS & NOZZLES



CO201-00-72	CO201-00-52	CO201-00-46	CO201-00-65	CO201-00-73
Long Femoral	Femoral	Femoral	Irrigation Femoral	Flexible
Spray Tip	Brush Tip	Spray Tip	Spray Tip	Tunnel Tip



CO201-00-80
Slash Shield



CO201-00-44	CO201-00-51
Fan Tip W/Splash	Fan
Shield	Spray Tip



CO201-00-37	CO201-00-50
Shower Spray Tip	Shower
W/Splash Shield	Spray Tip



ORDERING INFORMATION

KITS			
REF NUMBER	DESCRIPTION		
CO201-02-00	Pulse Lavage - Component Kit		
CO201-03-00	Pulse Lavage - Knee Kit (Shower Spray Tip with Splash Shield)		
CO201-04-00	Pulse Lavage - Hip Kit (Femoral Spray Tip & Spray Tip with Splash Shield)		
CO201-05-00	Pulse Lavage - Fan Spray (Fan Spray Tip with Splash Shield)		
TIPS			
CO201-00-37	Shower Spray Tip with Splash Shield		
CO201-00-44	Fan Spray Tip with Splash Shield		
CO201-00-46	Femoral Spray Tip		
CO201-00-50	Shower Spray Tip		
CO201-00-51	Fan Spray Tip		
CO201-00-52	Femoral Brush Tip		
CO201-00-65	Irrigation Femoral Spray Tip		
CO201-00-72	Long Femoral Spray Tip		
CO201-00-73	Flexible Tunnel Tip		
SPLASH SHIELD			
CO201-00-80	Splash Shield: 6 Inch (15.2cm)		

DISTRIBUTED BY



P: +61 426 246 267 E: cingularortho@gmail.com

www.cingular.com.au

CINGULAR ORTHOPAEDICS PTY LTD. All rights of distribution, also by photocopy, reprint (also in the form of extracts) or storing and recovery in any kind of data processing device are reserved and require our written approval. The information, illustrations, photography and general inaformation are correct and true at the time of print. However they cannot be regarded as binding warranty of quality. The user is responsible for checking the products prior to usewhether they are suited for the intended purpose. We can only be held liable for these instruments pursuant to our terms of delivery and payment. Subject to technical alterations. CING | PULSE LAVAGE: 2024 (V1)